

# U0008

## FUNCTIONS



Backrest adjustment : 0 - 68°



Thigh rest adjustment : 0 - 27°



Foot rest mechanical adjustment : 0 - 22°



Vascular position



Fowler position



Trendelenburg position : 0 - 12°

# HOSPITAL BED 2 MOTORS

## OPTIONS

### GROUP 1 / CASTOR OPTIONS



### GROUP 2 / SIDE RAIL OPTIONS



### GROUP 3 / ELECTRICAL ACCESSORY OPTIONS



### GROUP 4 / MATTRESS PLATFORM OPTIONS



### GROUP 5 / OTHER ACCESSORY OPTIONS



## DESCRIPTION OF OPTION CODES

<b>T16</b> CASTOR TWO WITH BRAKE (125 MM)	<b>E3</b> UNDER BED LIGHTING SYSTEM
<b>T4</b> METAL CASTOR TWO WITH BRAKE (125 MM)	<b>Y1</b> REMOVABLE 4 PARTED MESH PLATFORM.
<b>T14</b> CASTOR CENTRAL & DIRECTIONAL LOCK (125 MM)	<b>Y3</b> REMOVABLE DOUBLE WALLED PE BED PLATFORM
<b>T19</b> CASTOR CENTRAL & DIRECTIONAL LOCK (150 MM)	<b>W0</b> WITHOUT OTHER ACCESSORIES
<b>K0</b> WITHOUT SIDE RAIL	<b>W1</b> IV STAND : ONE HEIGHT ADJUSTMENT SYSTEM
<b>K4</b> WITH ALUMINIUM SIDE RAIL	<b>W2</b> LIFTING POLE : BEARING CAPACITY 75 Kg
<b>E0</b> WITHOUT ELECTRICAL ACCESSORIES	<b>W4</b> TELEVISION STAND FOR 17" LCD TV
<b>E1</b> BATTERY RECHARGEABLE LEAD GEL TYPE. OUTPUT VOLTAGE: 24 V DC. CAPACITY: 1,2 Ah.	<b>W5</b> TRACTION APPARATUS

**Must select only one of them.**

**Do not select the others if you select without accessory. If you select the others you can select all or 1-2-3 ... of them.**

## GENERAL SPECIFICATIONS

- Easy to clean four-piece bed platform.
- Head and foot movements with the help of an electrical motor.
- Trendelenburg movement by the help of shock absorber mechanism at foot side.
- Ergonomic, durable and easily cleanable aluminum rails system with maximum safety designed in accordance with the standards of EN 60601-2-52.
- Head and foot platforms with automatic sliding system.
- Easily removable and installable, hygienic, high resistant head board system designed of compact laminate and aluminum.
- Shock resistance plastic bumper wheels.
- Metal components dyed with anti-corrosion electrostatic powder paint after being treated with phosphate.
- Optional laminate color choice.<sup>(RLX)</sup>

## TECHNICAL SPECIFICATIONS

■ External dimension	: 102 x 222,5 x 49,5 cm
■ Width outside side rails	: 98 cm
■ Mattres platform	: 89 x 195 cm
■ Recommended mattress dimension	: 85 x 194 x 12 cm
■ Autoregression (Backrest / Legrest)	: 8 cm / 4 cm
■ Safe Operating Load (Inc. Mattres and Accessories)	: 230 kg
■ Weight	: ~ 101 kg
■ Volume	: 1,23 m <sup>3</sup>

## ELECTRICAL SPECIFICATIONS

■ Power Input	: 230 V AC 50 Hz Max 1 A
■ Power Output	: 24 V DC
■ Electrical Current Protection	: Class 1, Type B
■ Liquid and Dust Protection	: IP54 – IP66



\* RL : Laminate color options page 82.