

# U0011

# HOSPITAL BED WITH 2 CRANKS

## FUNCTIONS



Backrest adjustment : 0 - 70°



Thigh rest adjustment : 0 - 40°



Foot rest mechanical adjustment : 0 - 25°



Vascular position



Fowler position

## GENERAL SPECIFICATIONS

- Easily cleanable metal bed platform designed for the ventilation of the bed.
- Head and foot movements with rotating levers.
- Ergonomic, durable and easily cleanable metal rails system with maximum safety.
- Melamine coated chipboard head and foot panels.
- Metal components dyed with anti-corrosion electrostatic powder paint after being treated with phosphate.
- Optional melamine coated chipboard color choice.<sup>(RSX)</sup>

## TECHNICAL SPECIFICATIONS

- Width : 102 cm
- Length : 208 cm
- Height : 50 cm
- Weight : ~ 80 kg
- Volume : 0,64 m<sup>3</sup>

## OPTIONS

### GROUP 1 / CASTOR OPTIONS



### GROUP 2 / SIDE RAIL OPTIONS



### GROUP 3 / OTHER ACCESSORY OPTIONS



## DESCRIPTION OF OPTION CODES

- T0** WITHOUT CASTORS
- T4** METAL CASTOR TWO WITH BRAKE (125 MM)
- T16** CASTOR TWO WITH BRAKE (125 MM)
- K0** WITHOUT SIDE RAIL
- K1** WITH METAL SIDE RAIL
- K3** WITH METAL SIDE RAIL
- W0** WITHOUT OTHER ACCESSORIES
- W2** LIFTING POLE : BEARING CAPACITY 75 Kg
- W8** IV STAND : ONE HEIGHT ADJUSTMENT SYSTEM

Must select only one of them.

Do not select the others if you select without accessory.  
If you select the others you can select all or 1-2-3 ... of them.



\* RS : Melamine color options page 83.