

U0014

HOSPITAL BED WITH 1 CRANK

GENERAL SPECIFICATIONS

- Platform becomes from two sections. Bed floor is in the form of the grid so bed can receive the air.
- Head movement with rotating levers.
- Ergonomic, durable and easily cleanable metal rails system with maximum safety is optional.
- Portable head and foot boards.
- Shock resistance plastic bumper wheels.
- Metal components dyed with anti-corrosion electrostatic powder paint after being treated with phosphate.
- Optional melamine colour choice for the metal-framed wooden headboard.^(RSX)
- Optional laminate colour choice for the PE plastic headboard.^(RLX)

TECHNICAL SPECIFICATIONS


■ Width	: 105 cm
■ Length	: 204 cm
■ Height	: 50 cm
■ Backrest Adjustment	: 68°
■ Weight	: ~ 52 kg
■ Volume	: 0,64 m ³

OPTIONS

GROUP 1 / CASTOR OPTIONS



GROUP 2 / BOARD OPTIONS



GROUP 3 / SIDE RAIL OPTIONS



GROUP 4 / OTHER ACCESSORY OPTIONS



DESCRIPTION OF OPTION CODES

T0 WITHOUT CASTORS	W0 WITHOUT OTHER ACCESSORIES
T4 METAL CASTOR TWO WITH BRAKE (125 MM)	W2 LIFTING POLE : BEARING CAPACITY 75 Kg
T16 PLASTIC CASTOR TWO WITH BRAKE (125 MM)	W4 TELEVISION STAND FOR 17" LCD TV
B1 METAL BOARD	W5 TRACTION APPARATUS
B2 PE PLASTIC BOARD	W8 IV STAND : ONE HEIGHT ADJUSTMENT SYSTEM
K0 WITHOUT SIDE RAIL	■ Must select only one of them.
K1 WITH METAL SIDE RAIL	■ Do not select the others if you select without accessory.
K3 WITH METAL SIDE RAIL	■ If you select the others you can select all or 1-2-3 ... of them.



* RS : Melamine color options page 83.

* RL : Laminate color options page 82.