

U0017

BABY CODE

FUNCTIONS



Height : 85,5 cm



Trendelenburg position : 0 - 14°




Anti - trendelenburg position : 0 - 14°

OPTIONS

GROUP 1 / CASTOR OPTIONS



GROUP 2 / CABINET OPTIONS



DESCRIPTION OF OPTION CODES

- T0** WITHOUT CASTORS
- T7** CASTOR TWO WITH BRAKE (50 MM)
- D0** WITHOUT CABINET
- D1** METAL WITH CABINET

Must select only one of them.

GENERAL SPECIFICATIONS

- Manufacture from antibacterial, transparent flexi glass and portable.
- Trendelenburg and reverse trendelenburg position are adjusted by mechanical system.
- Epoxy painted steel construction is resistant to corrosion.
- Moves easily by 50 mm diameter four castors two of them with brake to fix the position. (Optional)

TECHNICAL SPECIFICATIONS

- | | | |
|------------------------------|---|---------------------|
| ■ Width | : | 45 cm |
| ■ Length | : | 79 cm |
| ■ Weight (Without Accessory) | : | ~ 11 kg |
| ■ Weight (Cabinet) | : | ~ 8 kg |
| ■ Volume | : | 0,34 m ³ |

