

U0043

## TROLLEY DRESSING



### GENERAL SPECIFICATIONS

- Easy cleanable and resistant ABS plastics top platform.
- At both sides there are one piece of quartet drug boxes produced from transparent plexiglass.
- Sliding shelf and drawer for drug preparation.
- There is serum containers made of ABS plastics.
- There is a stainless steel dustbin and dressing bowl.
- Height adjustable IV stand.
- Epoxy painted steel construction is, resistant to corrosion.
- Moves easily by 100 mm diameter four castors two of them with brake to fix the position.

### TECHNICAL SPECIFICATIONS

|          |                       |
|----------|-----------------------|
| ■ Width  | : 93 cm               |
| ■ Depth  | : 46 cm               |
| ■ Height | : 100 cm              |
| ■ Weight | : approx 47 kg        |
| ■ Volume | : 0,52 m <sup>3</sup> |

U0044

## TROLLEY DRUG DISTRIBUTION



### GENERAL SPECIFICATIONS

- Easy cleanable and resistant ABS plastics top platform.
- At both sides there are one piece of quartet and two each quinary drug boxes produced from transparent flexiglass.
- Sliding shelf and drawer for drug preparation.
- There is plastic dustbin one side of trolley and other side serum containers made of ABS plastic.
- Height adjustable IV stand.
- Epoxy painted steel construction is, resistant to corrosion.
- Moves easily by 100 mm diameter four castors two of them with brake to fix the position.

### TECHNICAL SPECIFICATIONS

|          |                       |
|----------|-----------------------|
| ■ Width  | : 93 cm               |
| ■ Depth  | : 46 cm               |
| ■ Height | : 100 cm              |
| ■ Weight | : approx 47 kg        |
| ■ Volume | : 0,52 m <sup>3</sup> |