

U0045

## TROLLEY EMERGENCY AND DRUG DISTRIBUTION

**GENERAL SPECIFICATIONS**

- Easy cleanable, strong ABS plastic upper platform and shelf for preparing medicine.
- Useful four pieces drawers can be locked.
- Stainless steel stand defibrillator.
- Stainless steel body.
- Long lived rail system.
- Easy moveable with 125 mm diameter casters and fixed each position with two casters having brake mechanisms.

**TECHNICAL SPECIFICATIONS**

■ Width	: 74,5 cm
■ Depth	: 63 cm
■ Height	: 105 cm
■ Weight	: approx 77 kg
■ Volume	: 0,61 m <sup>3</sup>

