

U0002 – HOSPITAL BED WITH 4 MOTORS

U0018 - MATTRESS



Images are for illustrative purposes only. The actual product may differ based on the configuration.

TECHNICAL SPECIFICATION

- Motorized beds for intensive care units and critical care units.
- Multifunctional type with 6-8 positions .
- Type: for adult patients.
- Type: with patient control. All buttons must be accompanied by graphic symbols.
- Functional features:
- Cardiac chair positioning.
- Trendelenburg positioning and Reverse-Trendelenburg positioning .
- Nurse control panel with patient control lock.
- Maximum load with label on the bed: ≥ 220 kg.
- Adjustable height: 78 - 79 cm.
- Wheels: Diameter: 12,5 cm. Lock system: all with seperate brake.
- Side guards consisting of 2 separate sections on each side of the bed (left, right).
- The rear section have a frame for an X-ray detector: $\geq 43 \times 43$ cm.
- Manual adjustment of important parts in case of power outage: Back Rest Motor With Manuel CPR Function.
- Power supply: 220 V, 50 Hz.
- Orthopedic mattresses with washable cover included.
- The mattress foam is 5 zone.